



HoTaq One-step Real-time RT-PCR Kit User Manual

384 Oyster Point Blvd, Suite 15
South San Francisco, CA 94080
Phone: 1 (888) MCLAB-88
Fax: 1 (650) 872-0253
www.mclab.com

Description

RT-PCR is widely used for measuring gene expression in tissue samples or cell culture systems. Traditionally, it is performed in separated two reaction steps. First-strand cDNA is reverse-transcribed from total RNA or poly (A)+ RNA using a reverse transcriptase. After then, the cDNA is amplified by PCR using a DNA polymerase in another reaction.

MCLAB the HoTaq One-step Real-time RT-PCR Kit offers a unique system for performing RT-PCR in a single step within a single tube with optimized reaction condition, which utilizes our own proprietary engineered highly purified thermostable McScript III reverse transcriptase and Taq DNA Polymerase. No additional reagents or steps are required once the reaction is initiated. This novel kit allows you to quantitatively detect specific RNA targets with high sensitivity, unparalleled convenience and wide dynamic range.

The technique reduces the risk of cross-contamination and minimizes reagents handling steps. This method is particularly useful for applications in which the expression of a small number of genes must be analyzed in many different total RNA samples, and robust amplifies high-abundance transcripts from crude total RNA preparations.

List of Components

Sufficient reagents are supplied in the MCLAB HoTaq One-step Real-time RT-PCR Kit for 200 rxns.

Upon receipt of the kit, immediately store the components at -20 °C in a freezer without a defrost cycle. It is recommended to prevent light exposure as little as possible.

2X HoTaq Real-time PCR Mix*: 1.25ml x 2

25X RT Enzyme Mix: 200ul

*HTRT400: Regular level of ROX, for Real-time PCR Machines ABI 7000, 7300, 7700, 7900.

HTRT400LR: Low level of ROX, for Real-time PCR Machines ABI 7500, Mx 3000P, Mx 3005P.

HTRT400RF: ROX Free, for Real-time PCR Machines BioRad iCycler MiniOpticon, Opticon 2, Chromo4, iQ5; Roche LightCycler 480; MJ Research DNA Engine Opticon 2, Chromo4; Corbett Roto-gene 3000, 6000.

Additional Materials Required

The following reagents, instruments and consumables are supplied by the user:

- Template RNA
- Gene-specific primers and probe
- DEPC-treated water
- Microcentrifuge
- Real-time thermal cycler
- PCR tubes/plates

One-step Real-time RT-PCR Procedure

1. Set up one-step real-time RT-PCR reaction in a PCR tube on ice as below:

Tube Cap Color	Component	Volume	Note
White	PCR Mix	12.5 µL	
Yellow	RT Mix	1 µL	
NA	Forward Primer	x	Final 50 to 900 nM*
NA	Reverse Primer	x	Final 50 to 900 nM*
NA	Probe	x	Final 50 to 250 nM*
NA	Template	x	10 pg to 1ug total RNA
NA	Nuclease-free water	To 25 µL	
Total Volume		25 µL	

* The concentration might need optimization for your specific targets.

For multiple reactions, master mix should be made with 5% extra reagents to reduce pipette error.

2. Gently mix thoroughly and then centrifuge briefly.

3. Place the tube or plate in the thermal cycler and set up program using associated software. Basic cycling conditions are as following:

One cycle at 50°C for 30 minutes;

One cycle at 95°C for 10 minutes;

Followed by 40* cycles of: 95°C for 15 seconds, 60°C * for 1 minute;

4°C hold (optional)

* cycle number and annealing temperature are experiment dependent.

Troubleshooting

Unexpected high or none Cq (quantification cycle) value

- Your RNA may be degraded. RNA integrity should be evaluated prior to the reaction.
- Your RNA sample might be contaminated with interferers for RT-PCR reaction. RNA purity should be evaluated prior to the reaction. Sample from crude total RNA preparation should be serial diluted as input and run multiple reactions. Re-precipitate or re-isolate RNA from the source if necessary.
- Poor primers and probe design. Be sure to follow guidelines when designing primers and probe for your target gene.
- Target gene may be at very weak expression level in your sample. Paralleled control gene amplification will be helpful for troubleshooting.

Reference

1. Higuchi R, et al (1992). BioTechnology 10:413
2. Higuchi R, et al (1993). BioTechnology 11:1026